



## PCR Grade Water

Cat. No.	Pack Size
water-00S	1,2 ml
water-25	25 ml
water-100	100 ml
water-500	500 ml

For *in vitro* use only

### Description:

PCR Grade Water is deionized and autoclaved water suitable for use in all experiments that require nuclease-free water. PCR Grade Water is prepared without chemical additives and it is Dnase, Rnase and nuclease-free.

### Application:

- All PCRs and RT-PCRs
- All other molecular biology techniques where pure water is needed

### Reagents Provided:

- PCR Grade Water

### Shipping and Storage conditions:

Shipping and storage at ambient temperature.

### Related products:

Product name	Pack size	Cat. No.
FIREPol® DNA Polymerase	500 U	01-01-00500
FIREPol® DNA Polymerase	1000 U	01-01-01000
FIREPol® DNA Polymerase	2000 U	01-01-02000
HOT FIREPol® DNA Polymerase	500 U	01-02-00500
HOT FIREPol® DNA Polymerase	1000 U	01-02-01000
dNTP MIX (20 mM of each)	20 µmol	02-31-00020
dNTP MIX (20 mM of each)	100 µmol	02-31-00100
dNTP SET (100 mM)	4 x 25 µmol	02-21-00100
dNTP SET (100 mM)	4 x 100 µmol	02-21-00400